

Cleaning Instructions for Humanscale Freedom Chair Textiles

All Freedom Textiles are treated with a stain resistant finish such as Teflon or Scotchguard. However, these finishes will not make the material completely impervious to stains or other invasive elements. Here is specific information for the removal of stains on the textile.

For Removal of Small Moist Stains: coffee, soda, juice, some salad dressings, etc.

- Combine 1 part detergent with 3 parts water to make a cleaning solution (1 Tbsp. of detergent with 3 Tbsp. of water.)
- Moisten soft sponge with solution
- Wet area with sponge and wipe it gently and thoroughly to remove stain

For Removal of Larger and Substance-based Stains: ketchup, mustard, salad dressing, etc.

- Remove excess substance with a knife edge or a coin
- Let stain dry (one or two days)
- Remove built-up stain with soft brush (a toothbrush works well)
- Using cleaning solution given above, wet area with sponge and wipe it gently and thoroughly to remove stain

Keeping Freedom Chair textiles clean and making them last:

- For liquid-based stains that penetrate the textile, removing the stain as quickly as possible will minimize the likelihood of permanent staining.
- For more substance-based stains that sit on top of the fabric, immediately scrape off the excess substance with a knife or coin.
- Always use water-soluble solvents such as dishwashing detergent. A few brand names are Ivory Liquid, Sun and Palmolive.
- Vacuuming the cushions/textiles on a weekly basis will extend the life of the material as well as identify any stains that need to be removed



available at:



800.531.3746

info@thehumansolution.com
thehumansolution.com